

# Oxylog® 2000 plus

Step up your performance with Oxylog® 2000 *plus*. The Oxylog® 2000 *plus* supports you in your daily challenge of saving peoples lives, no matter where the call takes you. Invasive or non-invasive, Oxylog® 2000 *plus* can meet this challenge by putting essential ventilation tools at your fingertips. The Oxylog® 2000 *plus* can make all the difference.



## TECHNICAL DATA

The Oxylog 2000 *plus* is a time-cycled, volume controlled emergency and transport ventilator with Pressure Support for patients requiring mandatory or assisted ventilation with a tidal volume from 100 mL upwards.

Dimensions (W x H x D)	285 x 184 x 175 mm / 11.10 x 7.24 x 6.89 inch (excluding handle)
Weight	Approximately 5.4 kg / 11.9 lbs (including internal battery)

### Gas supply

Supply gas	Medical Oxygen
Supply pressure	270 - 600 kPa at 100 L/min
Gas consumption for internal control	0.1 to 0.5 L/min

### Operating data

Ventilation Modes	VC-CMV, VC-AC, VC-SIMV, SpnCPAP
Options	Pressure Support Ventilation and Non-Invasive (mask) Ventilation
Special Functions	Apnea Ventilation (For switching over automatically to volume-controlled mandatory ventilation, if breathing stops)
Ventilation respiratory rate	2 to 50 /min ±1 /min (VC-SIMV) 5 to 50 /min ±1 /min (VC-CMV, VC-AC) 12 to 50 /min ±1 /min for apnea ventilation
Tidal volume Vt	100 to 2000 mL, BTPS*
Ventilation time ratio I:E (VC-CMV, VC-AC)	1:4 to 3:1
Inspiration time Ti (VC-SIMV, VC-SIMV / PS)	0.2 to 10 seconds
FiO <sub>2</sub> concentration	100% (No-AirMix) or approximately 40% (O <sub>2</sub> AirMix).**
PEEP	0 to 20 mbar / cmH <sub>2</sub> O
Trigger sensitivity (flow trigger)	3 to 15 L/min
Pressure support ΔP <sub>supp</sub>	0 to 35 mbar (relative to PEEP), slope adjustable in 3 steps
Maximum inspiratory flow	100 L/min (supply pressure > 350 kPa / 51 PSI, 80 L/min (supply pressure < 350 kPa / 51 PSI)
Measured value display	V <sub>Te</sub> , V <sub>e</sub> , M <sub>Vespon</sub> , RR, RR <sub>sp</sub> , PEEP, P <sub>mean</sub> , PIP, P <sub>plat</sub> , O <sub>2</sub>
Display	Technology: Electro-luminescence, Pixels: 240 * 128, Visible area: 108 * 56 mm / 4.25 x 2.20 inch



MT-4300-2007

**Oxylog® 2000 plus**  
Step up your performance

**Power supply**

Input voltage	19 V ±0.5 V DC
AC/DC power pack	Input: 100 to 240 V AC, Output: 19 V DC
DC/DC converter	Input: 12 / 24 / 28 V DC, Output: 19 V DC
Battery type	Lithium ion battery
Operating time (fully charged, "typical" ventilation)	Approximately 4 hours
Battery charging time	Approximately 5 hours

**Monitoring**

Supply pressure low	Supply pressure < 270 kPa / 39 PSI
Airway pressure (Paw high)	International: Adjustable from 20 to 60 mbar, USA: Adjustable from 20 to 100 cmH <sub>2</sub> O
Airway pressure (Paw low)	When pressure difference between inspiration and expiration < 5 mbar / cmH <sub>2</sub> O or when the set pressure level is not reached
Apnea alarm time Tapn	Adjustable from 15 to 60 seconds (not active when using NIV)
Leakage	VT <sub>e</sub> is approximately 40% lower than VT <sub>i</sub>
High respiratory rate	Patient breaths at a high spontaneous rate

**Operating Conditions**

Temperature	Temperature -20 to 50 °C / 14 to 122 °F
Atmospheric pressure	570 to 1200 hPa / 17 to 35 inches mmHg
Relative humidity	5 to 95%
Electromagnetic compatibility EMC	In accordance with ICE/EN 60601-1-2:2001 and ISO 10651-3
Airworthiness	In accordance with RTCA DO - 160D, sections 7, 8 & 21
Mechanical strength	In accordance with MIL STD 810F, method 514.5
Classification according to MDD 93/42/EEC	Class IIb
UMDNS-Code	18-098

\* BTPS: Body Temperature, Pressure, Saturated. Measured values referred to the conditions of the patient's lungs, body, temperature 37 °C / 99 °F, ambient pressure, water-vapor saturated gas.

\*\* Indirect measurement of O<sub>2</sub> concentration (calculated from two measured flows).

**HEADQUARTERS**

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

**REGION EUROPE CENTRAL  
AND EUROPE NORTH**

Dräger Medical GmbH  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

**REGION EUROPE SOUTH**

Dräger Médical S.A.S.  
Parc de Haute  
Technologie d'Antony 2  
25, rue Georges Besse  
92182 Antony Cedex, France  
Tel +33 1 46 11 56 00  
Fax +33 1 40 96 97 20  
[dlimfr-contact@draeger.com](mailto:dlimfr-contact@draeger.com)

**REGION MIDDLE EAST, AFRICA,  
CENTRAL AND SOUTH AMERICA**

Dräger Medical GmbH  
Branch Office Dubai  
Dubai Healthcare City, P.O. Box 505108  
Dubai, United Arab Emirates  
Tel + 971 436 24 762  
Fax + 971 436 24 761  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

**REGION ASIA / PACIFIC**

Dräger Medical  
South East Asia Pte Ltd  
25 International Business Park  
#04-27/29 German Centre  
Singapore 609916, Singapore  
Tel +65 6572 4388  
Fax +65 6572 4399  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

**REGION NORTH AMERICA**

Dräger Medical, Inc.  
3135 Quarry Road  
Telford, PA 18969-1042, USA  
Tel +1 215 721 5400  
Toll-free +1 800 437 2437  
Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

**Manufacturer:**

Dräger Medical Systems, Inc.  
Telford, PA 18969, USA  
The quality management system at  
Dräger Medical Systems, Inc. is  
certified according to ISO 13485,  
ISO 9001 and Annex II.3 of Directive  
93/42/EEC (Medical devices).